Adjusted Calcium (mg/dL)	
Calculation	Total Calcium + 0.8 x (4.0-Albumin)
Included In	Adjusted Calcium Phosphorus Product
	Comprehensive Metabolic Panel
	Comprehensive Metabolic Panel w/Phosphorus
Note	Calculation provided only when Albumin is <4.0 g/dL

Adjusted Calcium Phosphorus Product (mg²/dL²)	
Calculation	(Total Calcium (mg/dL) + 0.8 x (4.0-Albumin (g/dL)) x Phosphorus
Included In	Adjusted Calcium Phosphorus Product
	Comprehensive Metabolic Panel w/Phosphorus
Note	Calculation provided only when Albumin is <4.0 g/dL

A/G Ratio	
Calculation	Albumin/Globulin
Included In	Comprehensive Metabolic Panel
	Comprehensive Metabolic Panel w/Phosphorus

Anion Gap (mEq/L)	
Calculation	Sodium – (Chloride + CO2)
Included In	Basic Metabolic Panel
	Comprehensive Metabolic Panel
	Comprehensive Metabolic Panel w/ Phosphorus
	Electrolytes
	Renal Function Panel

Calcium Phosphorus Product (mg²/dL²)	
Calculation	Total Calcium x Phosphorus
Included In	Calcium Phosphorus Product
	Adjusted Calcium Phosphorus Product
	Comprehensive Metabolic Panel w/Phosphorus

eGFR (mL/min/1.73 m²)	
Calculation	175 x (Scr) ^{-1.154} x (age) ^{-0.203} x (0.742 if female) x (1.210 if African American)
Included In	eGFR (CKD, Non Dialysis)
Note	Estimated GFR (eGFR) using IDMS-Traceable Modification of Diet in Renal Disease (MDRD)

Globulin (g/dL)	
Calculation	Total Protein–Albumin
Included In	Comprehensive Metabolic Panel
	Comprehensive Metabolic Panel w/Phosphorus

Hemoglobin x 3 (g/dL)	
Calculation	Hemoglobin x 3
Included In	Complete Blood Count (CBC) & Differential
	Complete Blood Count (CBC) & Differential w/Reticulocytes
	Hemoglobin
	Hemoglobin & Hematocrit (H&H)
	Hemogram (Complete Blood Count w/o Differential)

Iron Status w/Iron & Transferrin	
Total Iron Binding Capacity (TIBC) (μg/dL)	
Calculation	Transferrin x 1.4
% Transferrin Saturation	
Calculation	(Iron/(Transferrin x 1.4)) x 100

Lipid Panel		
Cholesterol/HDL	Cholesterol/HDL Ratio	
Calculation	Cholesterol/HDL	
Included In	Lipid Panel (Coronary Risk Profile)	
Low Density Lipo	Low Density Lipoprotein (mg/dL)	
Calculation	Cholesterol – (Very Low Density Lipoprotein + HDL)	
Included In	Lipid Panel (Coronary Risk Profile)	
Very Low Density	Very Low Density Lipoprotein (mg/dL)	
Calculation	Triglycerides/5	
Included In	Lipid Panel (Coronary Risk Profile)	
Note	Only provided if Triglycerides are <400 mg/dL	

Prothrombin Time (Protime): INR	
Calculation	(PT Ratio) ^{ISI}
	PT Ratio = (Patient PT/Mean Normal PT) ^{ISI}
	Mean Normal PT = Geometric Mean
	ISI = International Sensitivity Index
Included In	Prothrombin Time (PT)

% Recirculation	
Calculation	(Systemic BUN – Arterial BUN)/(Systemic BUN – Venous BUN) x 100
Included In	Recirculation Study

HEMODIALYSIS CALCULATIONS	
Kt/V Equilibrated (eqKt/V)	
Calculation	(0.924 x lnKt/V) – ((0.395 x lnKt/V) / (Min / 60)) + 0.056
Included In	Kt/V Standard, URR
	Kt/V Standard, Natural Log, URR
Note	Leypoldt Formula
	For patient dialyzing 2 or 4-6 times per week
	For Kt/V Standard calculation purposes only, not reported

Kt/V Jindal	
Calculation	(0.04 x ((Pre BUN – Post BUN) / Pre BUN x 100) – 1.2)
Included In	Kt/V Jindal (Not K/DOQI Recommended)
Note	Jindal Formula
	The HD Adequacy Work Group feels this formula should not be used to measure delivered dose
	of Hemodialysis. (K/DOQI Clinical Practice Guidelines for Hemodialysis Adequacy: Update 2000,
	Guideline 2)

Kt/V Natural Log (InKt/V)	
Calculation	(-Ln((Post BUN/Pre BUN) - (0.008 x Treatment Time in mins/60)) + ((4-(3.5 x (Post BUN/Pre
	BUN)) x (Pre WT-Post WT)/Post WT))
Included In	Kt/V Natural Log, URR
	Kt/V Natural Log, URR, nPNA
	Kt/V Standard, Natural Log, URR
Note	Daugirdas II Formula
	The K/DOQI recommendations are: Prescribed dose of hemodialysis: Kt/V of 1.3
	Delivered dose of hemodialysis: Kt/V >1.2

Kt/V Standard (stdKt/V)	
Calculation	(168 x (1-exp(-eqKt/V)) / (Min/60)) / ((1-exp(-eqKt/V)) / eqKt/V + (168/Number of
	Treatment/(Min/60)) -1)
Included In	Kt/V Standard, URR
	Kt/V Standard, Natural Log, URR
Note	Leypoldt Formula
	For patient dialyzing 2 or 4-6 times per week

nPNA, Hemodialysis		
Calculation	1. Treatment #1: Beginning of week PNA (PCR) = Pre BUN /(36.3 + 5.48 x Kt/V Natural Log + 53.5/ Kt/V Natural Log) + 0.168	
	2. Treatment #2: Midweek PNA (PCR) = Pre BUN /(25.8 + 1.15 x Kt/V Natural Log + 56.4/ Kt/V Natural Log) + 0.168	
	3. Treatment #3: End of week PNA (PCR) = Pre BUN $/(16.3 + 4.3 \times \text{Kt/V Natural Log} + 56.6 / \text{Kt/V})$	
	Natural Log) + 0.168	
Included In	Kt/V Natural Log, URR, nPNA	
Note	nPNA calculation is only applicable to patients on thrice-weekly dialysis without significant	
	residual function. nPNA calculated from Kt/V without formal kinetic modeling according to	
	Depner T and Daugirdas J:JASN 1996:7:780-785.	

Urea Reduction Ratio (%)	
Calculation	(1 – (Post BUN/Pre BUN)) x 100
Included In	Kt/V Jindal (Not K/DOQI Recommended)
	Kt/V Natural Log, URR
	Kt/V Natural Log, URR, nPNA
	Kt/V Standard, Natural Log, URR
	Kt/V Standard, URR
	Urea Reduction Ratio w/Pre and Post BUN

Ultrafiltration Rate (UFR) (mL/kg/hr)	
Calculation	((pre-weight – post-weight) x 1000)/ (delivered time in mins/60)/ post-weight in kg
Included In	Kt/V Natural Log, URR
	Kt/V Natural Log, URR, nPNA
	Kt/V Standard, Natural Log, URR
	Kt/V Standard, URR

PD ADEQUACY	PD ADEQUACY CALCULATIONS	
Weekly Total Kt/V		
Calculation	Weekly Residual Kt/V + Weekly Dialysate Kt/V	
Weekly Residu	al Kt/V	
Calculation	((Urine Urea Nitrogen/BUN) x (Urine Volume (mL)/Urine Collection Time (min)) x 10.08) / VSA	
Note	Calculated if urine sample provided	
	· ·	
Weekly Dialysa	ate Kt/V	
Calculation	((Dialysate Urea Nitrogen/BUN) x (24 hour Dialysate Drain Volume (mL)/1000) x 7) / VSA	
Weekly Total C	Creatinine Clearance (Liters/week/1.73 m²)	
Calculation	Weekly Residual GFR + Weekly Dialysate Creatinine Clearance	
Weekly Residu	al GFR (Liters/week/1.73 m²)	
Calculation	Arithmetic Mean of Weekly Urea Clearance and Weekly Creatinine Clearance	
	((Urine Urea Nitrogen/BUN) x (Urine Volume (mL)/Urine Collection Time (min) x 10.08) + (Urine	
	Creatinine/Plasma Creatinine x Urine Volume (mL)/Urine Collection Time (min) x 10.08))/2 x	
	1.73/BSA	
Weekly Dialysa	ate Creatinine Clearance (Liters/week/1.73 m²)	
Calculation	(Dialysate Corrected Creatinine/Plasma Creatinine) x (24 hour Dialysate Drain Volume	
	(mL)/1000) x 7 x 1.73/BSA	
	tinine, 24 Hour (mg/dL)	
Calculation	Creatinine at 24 Hour Dwell – (Glucose at 24 Hour Dwell x 0.00010386)	
	eal Dialysis (g/kg/day)	
Calculation	(10.76 x ((0.69 x UNA) + 1.46)) / (VSA/0.58)	
	en Appearance (PNA) (g/day)	
Calculation	10.76 x ((0.69x UNA) + 1.46)	
UNA (g/day)		
Calculation	(24 Hour Drain Volume (mL) x 24 Hour Urea Dialysate)/100000 + (Urine Volume (mL) x Urine	
	Urea Nitrogen)/100000) x (1440/Total Urine Collection Time (min))	
Note	For PNA calculation purposes only, not reported	
Dadus C	Anna (DCA)	
Body Surface A	, · ·	
Calculation	Adult (≥16 years) uses DuBois and DuBois formula BSA (m²) = 0.007184 x Wt ^{0.425} x Ht ^{0.725}	
	Pediatric (< 16 years) uses Haycock formula BSA (m^2) = 0.024265 x Wt ^{0.5378} x Ht ^{0.3964}	
	where weight (Wt) is in kilograms and height (Ht) is in centimeters	
Valuma /\/\ f==	m Surface Area (Liters)	
	m Surface Area (Liters)	
Calculation	Adult (≥16 years) uses Hume and Weyers formula Malo: V= 14.012024 + 0.206785 × Wt + 0.104786 × Ht	
	Male: V= -14.012934 + 0.296785 x Wt + 0.194786 x Ht Female: V = -35.270121 + 0.183809 x Wt + 0.344547 x Ht	
	CIIIdle, V = -33.2/0121 + 0.103003 X W(+ 0.34434/ X T(

Pediatric (<16 years) uses Friis-Hansen formula

	$V = 0.135 \text{ x Wt}^{0.666} \text{ x Ht}^{0.535}$
	where weight (Wt) is in kilograms and height (Ht) is in centimeters
Included In	PD Adequacy
	Residual Urea Clearance (KrU)
Danita and Family	illhustion Tost (DET) Foot
	ilibration Test (PET) Fast
Calculation	tinine, 4 Hour (mg/dL) Creatinine at 4 Hour Dwell – (Glucose at 4 Hour Dwell x 0.00010386)
Calculation	Creatifilite at 4 hour Dwell – (Glucose at 4 hour Dwell x 0.00010386)
Corrected Crea	tinine D/P, 4 Hour
Calculation	Corrected Creatinine at 4 Hour Dwell/Plasma Creatinine at 2 Hour Dwell
Peritoneal Equi	ilibration Test (PET) Standard
Corrected Crea	tinine, 0 Hour, 2 Hour, 4 Hour (mg/dL)
Calculation	Creatinine at 0 or 2 or 4 Hour Dwell – (Glucose at 0 or 2 or 4 Hour Dwell x 0.00010386)
<u> </u>	7
	tinine D/P, 0 Hour
Calculation	Corrected Creatinine at 0 Hour Dwell/Plasma Creatinine at 2 Hour Dwell
Corrected Crea	tinine D/P, 2 Hour
Calculation	Corrected Creatinine at 2 Hour Dwell/Plasma Creatinine at 2 Hour Dwell
Calculation	Corrected Creatiffine at 2 flour Dwell/Flasifia Creatiffine at 2 flour Dwell
Corrected Crea	tinine D/P, 4 Hour
Calculation	Corrected Creatinine at 4 Hour Dwell/Plasma Creatinine at 2 Hour Dwell
Glucose D/D0,	2 Hour
Calculation	Glucose at 2 Hour Dwell/Glucose at 0 Hour Dwell
Glucose D/D0,	
Calculation	Glucose at 4 Hour Dwell/Glucose at 0 Hour Dwell
Urea D/P, 0 Ho	ur .
Calculation	Urea at 0 Hour Dwell/Urea at 2 Hour Dwell
Carcaration	orea at o riodi bwelly orea at 2 riodi bwell
Urea D/P, 2 Ho	ur
Calculation	Urea at 2 Hour Dwell/Urea at 2 Hour Dwell
Urea D/P, 4 Ho	ur
Calculation	Urea at 4 Hour Dwell/Urea at 2 Hour Dwell
Fluid, 24-Hour	
Calculation	tinine, 24 Hour (mg/dL) Creatinine at 24 Hour Dwell – (Glucose at 24 Hour Dwell x 0.00010386)
CalculatiOII	Creatifine at 24 flour Dwell – (Glucose at 24 flour Dwell x 0.00010300)
Fluid, Overnigh	t Dwell
	tinine, Overnight (mg/dL)
Calculation	Creatinine Overnight Dwell – (Glucose Overnight Dwell x 0.00010386)
	· · · · · · · · · · · · · · · · · · ·

Creatinine Clearance, Urine (Residual Renal Creatinine Clearance)	
Calculation	(Urine Creatinine/Blood Creatinine) x (Urine Volume (mL)/Urine Collection Time (min)) x (1.73/BSA)
Included In	Urine Creatinine Clearance

Body Surface Area (BSA)	
Calculation	Adult (≥16 years) uses DuBois and DuBois formula
	BSA (m^2) = 0.007184 x Wt ^{0.425} x Ht ^{0.725}
	Pediatric (< 16 years) uses Haycock formula
	BSA (m^2) = 0.024265 x Wt ^{0.5378} x Ht ^{0.3964}
	where weight (Wt) is in kilograms and height (Ht) is in centimeters

Residual Urea Clearance, KrU – for Hemodialysis only	
KrU (mL/min)	
Calculation	(Urine Urea Nitrogen x Urine Volume (mL)) / (Blood BUN x 0.9 x Total Urine Collection Time
	(min))

Kt/V Residual		
Calculation	(Urine Urea Nitrogen/Blood BUN) x (Urine Volume (mL)/ Urine Collection Time (min)) x	
	(10.08/VSA)	

Volume from Surface Area (Liters)		
Calculation	Adult (≥16 years) uses Hume and Weyers formula	
	Male: V= -14.012934 + 0.296785 x Wt + 0.194786 x Ht	
	Female: V = -35.270121 + 0.183809 x Wt + 0.344547 x Ht	
	Pediatric (<16 years) uses Friis-Hansen formula	
	$V = 0.135 \times Wt^{0.666}0.666 \times Ht^{0.535}$	
	where weight (Wt) is in kilograms and height (Ht) is in centimeters	
Included In	PD Adequacy	
	Residual Urea Clearance (KrU)	

Urine Creatinine, 24 Hour (mg/24 hour)		
Calculation	(Urine Creatinine/100) x Total Urine Volume in mL	
Included In	24 Hour Urine Creatinine	

Urine Urea Nitrogen, 24 Hour (g/24 hour)		
Calculation	((Urine Urea Nitrogen in mg/dL/100) x Total Urine Volume in mL) / 1000	
Included In	24 Hour Urine Urea Nitrogen	