

Adjusted Calcium (mg/dL)	Total Calcium + 0.8 x (4.0-Albumin)
Included In	Adjusted Calcium Phosphorus Product
	Comprehensive Metabolic Panel
	Comprehensive Metabolic Panel w/Phosphorus
Note	Calculation provided only when Albumin is <4.0 g/dL

Adjusted Calcium Phosphorus	(Total Calcium (mg/dL) + 0.8 x (4.0-Albumin (g/dL)) x Phosphorus
Product (mg ² /dL ²)	
Included In	Adjusted Calcium Phosphorus Product
	Comprehensive Metabolic Panel w/Phosphorus
Note	Calculation provided only when Albumin is <4.0 g/dL

A/G Ratio	Albumin/Globulin
Included In	Comprehensive Metabolic Panel
	Comprehensive Metabolic Panel w/Phosphorus

Anion Gap (mEq/L)	Sodium – (Chloride + CO2)
Included In	Basic Metabolic Panel
	Comprehensive Metabolic Panel
	Comprehensive Metabolic Panel w/ Phosphorus
	Electrolytes
	Renal Function Panel

Calcium Phosphorus Product	Total Calcium x Phosphorus
(mg^2/dL^2)	
Included In	Calcium Phosphorus Product
	Adjusted Calcium Phosphorus Product
	Comprehensive Metabolic Panel w/Phosphorus

eGFR (mL/min/1.73 m²)	175 x (Scr) ^{-1.154} x (age) ^{-0.203} x (0.742 if female) x (1.210 if African American)
Included In	eGFR (CKD, Non Dialysis)
Note	Estimated GFR (eGFR) using IDMS-Traceable Modification of Diet in Renal Disease (MDRD)

Globulin (g/dL)	Total Protein–Albumin
Included In	Comprehensive Metabolic Panel
	Comprehensive Metabolic Panel w/Phosphorus

Hemoglobin x 3 (g/dL)	Hemoglobin x 3
Included In	Complete Blood Count (CBC) & Differential
	Complete Blood Count (CBC) & Differential w/Reticulocytes
	Hemoglobin
	Hemoglobin & Hematocrit (H&H)
	Hemogram (Complete Blood Count w/o Differential)

Page 1 of 7 Rev. January 2022



Transferrin x 1.4
(Iron/(Transferrin x 1.4)) x 100
Cholesterol/HDL
Cholesterol – (Very Low Density Lipoprotein + HDL)
Triglycerides/5
Only provided if Triglyceride is <400 mg/dL
(PT Ratio) ^{ISI}
PT Ratio = (Patient PT/Mean Normal PT) ^{ISI}
Mean Normal PT = Geometric Mean
ISI = International Sensitivity Index
(Systemic BUN – Arterial BUN)/(Systemic BUN – Venous BUN) x 100

Page 2 of 7 Rev. January 2022



HEMODIALYSIS CALCULATIONS	
Kt/V Equilibrated (eqKt/V)	(0.924 x lnKt/V) – ((0.395 x lnKt/V) / (Min / 60)) + 0.056
Included In	Kt/V Standard, URR
	Kt/V Standard, Natural Log, URR
Note	Leypoldt Formula
	For patient dialyzing 2 or 4-6 times per week
	For Kt/V Standard calculation purposes only, not reported
Kt/V Jindal	(0.04 x ((Pre BUN – Post BUN) / Pre BUN x 100) – 1.2)
Included In	Kt/V Jindal (Not K/DOQI Recommended)
Note	Jindal Formula
	The HD Adequacy Work Group feels this formula should not be used to measure
	delivered dose of Hemodialysis. (K/DOQI Clinical Practice Guidelines for Hemodialysis
	Adequacy: Update 2000, Guideline 2)
Kt/V Natural Log (InKt/V)	(-Ln((Post BUN/Pre BUN) - (0.008 x Treatment Time in mins/60)) + ((4-(3.5 x (Post
Kty v Ivatural Log (IIIKty v)	BUN/Pre BUN)) x (Pre WT-Post WT)/Post WT))
Included In	Kt/V Natural Log, URR
	Kt/V Natural Log, URR, nPNA
	Kt/V Standard, Natural Log, URR
Note	Daugirdas II Formula
	The K/DOQI recommendations are: Prescribed dose of hemodialysis: Kt/V of 1.3
	Delivered dose of hemodialysis: Kt/V >1.2
Kt/V Standard (stdKt/V)	(168 x (1-exp(-eqKt/V)) / (Min/60)) / ((1-exp(-eqKt/V)) / eqKt/V + (168/Number of
ito, i cominanta (comito, i)	Treatment/(Min/60)) -1)
Included In	Kt/V Standard, URR
	Kt/V Standard, Natural Log, URR
Note	Leypoldt Formula
	For patient dialyzing 2 or 4-6 times per week
nPNA, Hemodialysis	1. Treatment #1: Beginning of week PNA (PCR) = Pre BUN /(36.3 + 5.48 x Kt/V Natural
,	Log + 53.5/ Kt/V Natural Log) + 0.168
	2. Treatment #2: Midweek PNA (PCR) = Pre BUN /(25.8 + 1.15 x Kt/V Natural Log +
	56.4/ Kt/V Natural Log) + 0.168
	3. Treatment #3: End of week PNA (PCR) = Pre BUN /(16.3 + 4.3 x Kt/V Natural Log +
	56.6/ Kt/V Natural Log) + 0.168
Included In	Kt/V Natural Log, URR, nPNA
Note	nPNA calculation is only applicable to patients on thrice-weekly dialysis without
	significant residual function. nPNA calculated from Kt/V without formal kinetic modeling
	according to Depner T and Daugirdas J:JASN 1996:7:780-785.
Urea Reduction Ratio (%)	(1 – (Post BUN/Pre BUN)) x 100
Included In	Kt/V Jindal (Not K/DOQI Recommended)
	Kt/V Natural Log, URR
	Kt/V Natural Log, URR, nPNA
	Kt/V Standard, Natural Log, URR
	Kt/V Standard, URR
	Urea Reduction Ratio w/Pre and Post BUN
Ultrafiltration Rate (UFR)	((pre-weight – post-weight) x 1000)/ (delivered time in mins/60)/ post-weight in kg
(mL/kg/hr)	
Included In	Kt/V Natural Log, URR
	Kt/V Natural Log, URR, nPNA
	Kt/V Standard, Natural Log, URR
	Kt/V Standard, URR

Page **3** of **7** Rev. January 2022



PD ADEQUACY CALCULATIONS	
Weekly Total Kt/V	Weekly Residual Kt/V + Weekly Dialysate Kt/V
Weekly Total Kt/ V	Weekly Residual KL/ V + Weekly Dialysate KL/ V
Weekly Residual Kt/V	((Urine Urea Nitrogen/BUN) x (Urine Volume (mL)/Urine Collection Time (min)) x 10.08) / VSA
Note	Calculated if urine sample provided
Weekly Dialysate Kt/V	((Dialysate Urea Nitrogen/BUN) x (24 hour Dialysate Drain Volume (mL)/1000) x 7) / VSA
Weekly Total CrCl (Liters/week/1.73 m²)	Weekly Residual GFR + Weekly Dialysate Creatinine Clearance
Weekly Residual GFR (Liters/week/1.73 m²)	Arithmetic Mean of Weekly Urea Clearance and Weekly Creatinine Clearance ((Urine Urea Nitrogen/BUN) x (Urine Volume (mL)/Urine Collection Time (min) x 10.08) + (Urine Creatinine/Plasma Creatinine x Urine Volume (mL)/Urine Collection Time (min) x 10.08))/2 x 1.73/BSA
Weekly Dialysate CrCl (Liters/week/1.73 m²)	(Dialysate Corrected Creatinine/Plasma Creatinine) x (24 hour Dialysate Drain Volume (mL)/1000) x 7 x 1.73/BSA
Weekly Residual CrCl (Liters/week/1.73 m²)	(Urine Creatinine/Blood Creatinine) x (Urine Volume (mL)/Urine Collection Time (min)) x (1.73/BSA) x 10.08
Creatinine Clearance (mL/min/1.73m2)	(Urine Creatinine/Blood Creatinine) x (Urine Volume (mL)/Urine Collection Time (min)) x (1.73/BSA)
Corrected Creatinine, 24 Hour (mg/dL)	Creatinine at 24 Hour Dwell – (Glucose at 24 Hour Dwell x 0.00010386)
nPNA, Peritoneal Dialysis (g/kg/day)	(10.76 x ((0.69 x UNA) + 1.46)) / (VSA/0.58)
Protein Nitrogen Appearance (PNA) (g/day)	10.76 x ((0.69x UNA) + 1.46)
UNA (g/day)	(24 Hour Drain Volume (mL) x 24 Hour Urea Dialysate)/100000 + (Urine Volume (mL) x Urine Urea Nitrogen)/100000) x (1440/Total Urine Collection Time (min))
Note	For PNA calculation purposes only, not reported
Body Surface Area (BSA)	Adult (≥16 years) uses DuBois and DuBois formula BSA (m²) = 0.007184 x Wt ^{0.425} x Ht ^{0.725} Pediatric (< 16 years) uses Haycock formula BSA (m²) = 0.024265 x Wt ^{0.5378} x Ht ^{0.3964} where weight (Wt) is in kilograms and height (Ht) is in centimeters
Volume (V) from Surface Area (Liters)	Adult (≥16 years) uses Hume and Weyers formula Male: V= -14.012934 + 0.296785 x Wt + 0.194786 x Ht Female: V = -35.270121 + 0.183809 x Wt + 0.344547 x Ht Pediatric (<16 years) uses Friis-Hansen formula V = 0.135 x Wt ^{0.666} x Ht ^{0.535} where weight (Wt) is in kilograms and height (Ht) is in centimeters

Page 4 of 7 Rev. January 2022



Peritoneal Equilibration Test (PET)	Peritoneal Equilibration Test (PET) Fast		
Corrected Creatinine, 4 Hour (mg/dL)	Creatinine at 4 Hour Dwell – (Glucose at 4 Hour Dwell x 0.00010386)		
Corrected Creatinine D/P, 4 Hour	Corrected Creatinine at 4 Hour Dwell/Plasma Creatinine		
Peritoneal Equilibration Test (PET)	Standard		
Corrected Creatinine, 0 Hour, 2 Hour, 4 Hour (mg/dL)	Creatinine at 0 or 2 or 4 Hour Dwell – (Glucose at 0 or 2 or 4 Hour Dwell x 0.00010386)		
Corrected Creatinine D/P, 0 Hour	Corrected Creatinine at 0 Hour Dwell/Plasma Creatinine		
Corrected Creatinine D/P, 2 Hour	Corrected Creatinine at 2 Hour Dwell/Plasma Creatinine		
Corrected Creatinine D/P, 4 Hour	Corrected Creatinine at 4 Hour Dwell/Plasma Creatinine		
Glucose D/D0, 2 Hour	Glucose at 2 Hour Dwell/Glucose at 0 Hour Dwell		
Glucose D/D0, 4 Hour	Glucose at 4 Hour Dwell/Glucose at 0 Hour Dwell		
Urea D/P, 0 Hour	Urea at 0 Hour Dwell/Plasma Urea		
Urea D/P, 2 Hour	Urea at 2 Hour Dwell/Plasma Urea		
Urea D/P, 4 Hour	Urea at 4 Hour Dwell/Plasma Urea		
Peritoneal Equilibration Test (PET)	Modified		
Corrected Creatinine, 0 Hour, 1 Hour, 2 Hour, 4 Hour (mg/dL)	Creatinine at 0 or 1 or 2 or 4 Hour Dwell – (Glucose at 0 or 1 or 2 or 4 Hour Dwell x 0.00010386)		
Corrected Creatinine D/P, 4 Hour	Corrected Creatinine at 4 Hour Dwell/Plasma Creatinine		
Sodium D/P, 0 Hour	Sodium at 0 Hour Dwell / Plasma Sodium		
Sodium D/P, 1 Hour	Sodium at 1 Hour Dwell / Plasma Sodium		
Sodium D/P, 2 Hour	Sodium at 2 Hour Dwell / Plasma Sodium		
Sodium D/P, 4 Hour	Sodium at 4 Hour Dwell / Plasma Sodium		
Fluid, 24-Hour Dwell			
Corrected Creatinine, 24 Hour (mg/dL)	Creatinine at 24 Hour Dwell – (Glucose at 24 Hour Dwell x 0.00010386)		
[
Fluid, Overnight Dwell Corrected Creatinine, Overnight (mg/dL)	Creatinine Overnight Dwell – (Glucose Overnight Dwell x 0.00010386)		
· · ·			

Page 5 of 7 Rev. January 2022



CALCOLATIONS	
24 Hour Urine Creatinine Clearance (Residual Renal Creatinine Clearance)	
Creatinine Clearance	(Urine Creatinine/Blood Creatinine) x (Urine Volume (mL)/Urine Collection Time (min)) x
(mL/min/1.73m ²)	(1.73/BSA)
Body Surface Area (BSA)	Adult (≥16 years) uses DuBois and DuBois formula
	BSA (m ²) = $0.007184 \times Wt^{0.425} \times Ht^{0.725}$
	Pediatric (< 16 years) uses Haycock formula
	BSA (m ²) = $0.024265 \times Wt^{0.5378} \times Ht^{0.3964}$
	where weight (Wt) is in kilograms and height (Ht) is in centimeters
	where weight (wt) is in knograms and height (itt) is in centimeters
Residual Urea Clearance, KrU – fo	ar Hemodialysis only
KrU (mL/min)	(Urine Urea Nitrogen x Urine Volume (mL)) / (Blood BUN x 0.9 x Total Urine Collection
KIO (IIIL/IIIIII)	, ,,,,
	Time (min))
when it is	
Kt/V Residual	(Urine Urea Nitrogen/Blood BUN) x (Urine Volume (mL)/ Urine Collection Time (min)) x
	(10.08/VSA)
6.3.4	T
Volume (V) from Surface Area	Adult (≥16 years) uses Hume and Weyers formula
(Liters)	Male: V= -14.012934 + 0.296785 x Wt + 0.194786 x Ht
	Female: V = -35.270121 + 0.183809 x Wt + 0.344547 x Ht
	Pediatric (<16 years) uses Friis-Hansen formula
	$V = 0.135 \times Wt^{0.666}0.666 \times Ht^{0.535}$
	where weight (Wt) is in kilograms and height (Ht) is in centimeters
24 Hour Urine Creatinine	
Urine Creatinine, 24 Hour	((Urine Creatinine in mg/dL x Urine Volume in mL)/100) x (1440 / Total Urine Collection
(mg/24 hr)	Time in min)
24 Hour Urine Urea Nitrogen	
Urine Urea Nitrogen, 24 Hour	((Urine Urea Nitrogen in mg/dL/100) x (Urine Volume in mL/(Total Urine Collection Time
(g/24 hr)	in mins/1440)))/1000
,	
24 Hour Urine Total Protein with	Creatinine
Urine Protein, 24 Hour	(Urine Total Protein/100) x ((Urine Volume in mL)/(Total Urine Collection Time in
(mg/24 hr)	mins/1440))
(6) =)	
Urine Total Protein/Creatinine	((Urine Total Protein/100) x ((Urine Volume in mL)/(Total Urine Collection Time in
Ratio, 24 Hour	mins/1440)))/(((Urine Creatinine/100) x (Urine Volume in mL/(Total Urine Collection
(mg/g creat)	Time in mins /1440)))/1000)
(iiig/g cicat/	Time in timis / 1440//// 1000/
Uring Creatining 24 Usus	((Urine Creatinine/100) x (Urine Volume in mL/(Total Urine Collection Time in
Urine Creatinine, 24 Hour	
(g/24 hr)	mins/1440)))/1000
Random Urine Total Protein with	
Urine Total Protein/Creatinine	(Total Protein, Random Urine/Creatinine, Random Urine) x 1000
Ratio (mg/g creat)	

Page **6** of **7** Rev. January 2022



Gotch PD QA		
Corrected Creatinine, PD QA (mg/dL)	(Corrected Creatinine, PD QA – Glucose PD QA x 0.00010386)	
Total Protein, PD QA (g/dL)	Total Protein, PD QA / 1000	
Gotch PD Exchange 1		
Corrected Creatinine, PD Exchange	(Corrected Creatinine, PD Exchange 1 – Glucose PD Exchange 1 x 0.00010386)	
1 (mg/dL)	(corrected diseasement) is Excitating 1 x disease is Excitating 1 x disease)	
Total Protein, PD Exchange 1	Total Protein, PD Exchange 1 / 1000	
(g/dL)	January Committee Committe	
Gotch PD Exchange 2		
Corrected Creatinine, PD Exchange 2 (mg/dL)	(Corrected Creatinine, PD Exchange 2 – Glucose PD Exchange 2 x 0.00010386)	
Total Protein, PD Exchange 2	Total Protein, PD Exchange 2 / 1000	
(g/dL)	Total Frotein, FB Exchange 27 1000	
Gotch PD Exchange 3		
Corrected Creatinine, PD Exchange 3 (mg/dL)	(Corrected Creatinine, PD Exchange 3 – Glucose PD Exchange 3 x 0.00010386)	
Total Protein, PD Exchange 3 (g/dL)	Total Protein, PD Exchange 3 / 1000	
Gotch PD Exchange 4		
Corrected Creatinine, PD Exchange 4 (mg/dL)	(Corrected Creatinine, PD Exchange 4 – Glucose PD Exchange 4 x 0.00010386)	
Total Protein, PD Exchange 4	Total Protein, PD Exchange 4 / 1000	
(g/dL)	Total Frotein, FD Exchange 47 1000	
Gotch PD Exchange 5		
Corrected Creatinine, PD Exchange 5 (mg/dL)	(Corrected Creatinine, PD Exchange 5 – Glucose PD Exchange 5 x 0.00010386)	
<u> </u>		
Total Protein, PD Exchange 5 (g/dL)	Total Protein, PD Exchange 5 / 1000	

Page **7** of **7** Rev. January 2022