Specimen Integrity | Avoid Canceled Tests

Label & Collection

CONFIRM correct patient correct collection date SMITH, JOHN 07/07/21 Use fill line/range as a visual guide **CONFIRM** correct tube type

Centrifuge Gel Tubes

Process all specimens within 2 hours.

- Centrifuge before refrigerating.
- Red Top tube requires minimum ½ hr to clot after draw.
- Opposing tube carriers must be identical.
- Centrifuge full cycle.
- Inspect tubes for proper separation.




Properly Spun Tube

Plasma/Serum (clear yellow part) contains no red blood cells.

Gel is firmly wedged against sides, separating plasma from red blood cells.

Red blood cells are below gel barrier.



Common Errors

Unspun Tube

- Gel on bottom.
- No separation of red blood cells from plasma/serum.

Spin before shipping



 Separation of red blood cells from plasma/serum.

> Spin immediately

Poorly Spun

- Plasma/serum will be pink or red.
- No complete gel barrier.

Properly balance centrifuge



Clotted

Clots adhere to tube when tilted.

Invert tube **IMMEDIATELY** after draw



requirement not met. Use nominal fill range/line as a

visual guide

Hemolysis • Plasma/serum will be

pink or red.

Follow proper specimen collection

Properly Packaged Specimens

Refrigerated

Shipment does not contain needles.


Ice Packs placed properly, using correct size/quantity.

- Small/Medium Box = 2
- Large Box = 4

Refrigerated Specimens packaged in biohazard transport bag with absorbent sheet and does not contain room temperature specimens.

Foam Layers placed on top and bottom of bagged specimens to protect from direct contact with ice packs.

Room Temperature



Shipment does not contain needles or ice packs.

Room Temperature Specimens packaged in biohazard transport bag with absorbent sheet and does not contain refrigerated specimens.

Foam Layers placed on top and bottom of bagged specimens.



Client Services **800.800.5655** Visit www.ascend-u.com for further training.