





Environmental Samples


Collection & Handling

Proper sterile collection technique is required when collecting environmental samples.

	Type	Tests	Collection	Refrigerate
	Sterile Cup w/Transfer Port	<i>For collection only</i>	<ul style="list-style-type: none"> Discard after collection Do not ship 	—
	Yellow Stopper	<ul style="list-style-type: none"> LAL Colony Count Ultrapure LAL 	<ul style="list-style-type: none"> Collect sample using Sterile Cup w/Transfer Port Fill tube until flow ceases 	✓
	Conical Tube	<ul style="list-style-type: none"> Water Analysis Conductivity 	<ul style="list-style-type: none"> Fill to MAX FREEZE line Secure lid 	✓
	Black Top	<ul style="list-style-type: none"> Electrolyte Dialysate Plus Glucose, Dialysate 	<ul style="list-style-type: none"> Collect sample using Sterile Cup w/Transfer Port Fill tube until flow ceases 	✓
	Sterile Cup	<ul style="list-style-type: none"> UltraPure Colony Count 	<ul style="list-style-type: none"> Collect at least 100 mL Secure lid 	✓

Labeling & Packaging

Verify label



CONFIRM correct source

CONFIRM correct collection date

CONFIRM correct tube type

MACHINE 03, DIALYSATE
 Mm: ENV80969 Dob: 01/01/2020 Fac: ABC
 06/01/21
 CT1
 12345678910
 YELLOW STOPPER

