






Order of Draw	Mix by Inverting	Centrifuge	Refrigerate
 ROYAL BLUE	8-10 times	No	Yes
 LIGHT BLUE	4 times	No	No
 RED	5-10 times	Yes (after clots)	Yes
 GREEN	8-10 times	Yes	Yes
 TALL GREEN	8-10 times	Yes	Yes
 LAVENDER	8-10 times	No	Yes

#### POST ORDER OF DRAW



GOLD/WHITE

8-10 times

Yes

Yes



Example of  
one inversion

Draw all specimens, except **POST BUN**  
(Gold/White), prior to administering heparin