

CONFIRM correct patient

CONFIRM correct collection date

CONFIRM correct tube

Process All Specimens within 2 hours

**Follow the
Order of Draw**
**Immediately
Mix by Inverting**
**Run Balanced
Centrifuge**
Refrigerate

ROYAL BLUE

8-10 times


LIGHT BLUE

4 times


RED/YELLOW

5-10 times


After clotting

GREEN

8-10 times


Immediately

LAVENDER

8-10 times


GOLD/WHITE

8-10 times


Immediately

**GREEN/YELLOW
ALTERNATE**
POST